



1106 Market Street
Pocomoke City, MD 21851
sales@newtownehyperbarics.com
410-575-4220

PRODUCT DATA SHEET– Newtowne Shoe

For over 20 years, Newtowne Hyperbarics has been a leader in designing and manufacturing mild hyperbaric chambers for homes. Our 100% American-built machines are suitable for home or clinical use. Our Class 4 Chamber has stood the test of time, with thousands being used in medical offices and homes all over America. A 2 Year Warranty is included with the ability to add 3 additional years of full warranty for every chamber new or used.

At Newtowne Hyperbarics all of our chambers are designed and manufactured with your safety as our priority. We can boast a 100% safety record since we began manufacturing chambers in 2001. All of our chambers come equipped with the following safety features:

- ◆ Dual, redundant pressure relief valves : only one is required for proper ventilation, the second is installed for safety
- ◆ Safety override relief valve on the compressor: prevents chamber from over inflating.
- ◆ Medical grade materials including the chamber material and components.

**Our products are cleared by and registered with the FDA
FDA 510k #K051759 - FDA OWNER OPERATOR #10037097**

**Sale or use of this device requires a prescription.
Federal law restricts this device for sale only by or on the order of a
physician, dentist, or other qualified professional.**

Frequency and Duration of Chamber Sessions
Chamber treatment schedules are established by the physician , dentist or therapist.

Indications for use

Mild Hyperbaric Chambers are intended for the treatment of Acute Mountain Sickness. All other conditions are considered "off label".

Environments for use

Chambers may be used as normal in inside or outdoor (-25° F to 120° F) environs. Chamber should operate as normal regardless of humidity or altitude.

PRODUCT DATA SHEET—NEWTOWNE SHOE



- ◆ 4psi / 1.3 ATA operation
- ◆ 80L Compressor
- ◆ 2 viewports
- ◆ 2 extra valve ports for accessories
- ◆ Shipping included / no tax
- ◆ Internal gauge \$150 upgrade
- ◆ No chair or matt provided: a legless gaming chair or 6 foot mat is recommended.

STANDARD SHOE DIMENSIONS

Fully Inflated:

68" long / 37" wide / 42" tall

FRAME DIMENSIONS:

64.5" long / 27.75" wide / 42" high—Head End/
24" high—Foot End

LONG SHOE DIMENSIONS

Fully Inflated:

81" long / 37" wide / 43" tall

FRAME DIMENSIONS:

68" long / 27.75" wide / 43" high—
Head End/ 24" high—Foot End

Chamber foot can be dropped down for easier access whether using on the floor or on a platform/bed



PRODUCT DATA SHEET—NEWTOWNE SHOE

CHAIRS

The Newtowne Shoe does not come with a mat or chair. We recommend using this chamber with a legless type of chair or mattress. Dimensions should be 24-25" wide, the length can be up to 64" in the shoe and 75" in the Long Shoe.

Here are some we have tried at the factory that can be purchased on AMAZON:

